

## Mycoplasma-ExS Efficacy Summary

Activity Claim	Test Protocol	Contact Time	Test Organisms
Bactericidal	EN13697	5 minutes	Pseudomonas aeruginosa, Escherichia coli, Enterococcus hirae, Staphylococcus aureus

## Spore-EX Efficacy Summary

Activity Claim	Test Protocol	Contact Time	Test Organisms
Bactericidal	EN1276	1 minute	Escherichia coli, Staphylococcus aureus, Pseudomonas aeruginosa, Salmonella typhimurium, Enterobacter cloaca, Enterococcus hirae, Lactobacillus planetarium, Klebsiella pneumonia, Enterococcus faecalis, Proteus mirabilis, Streptococcus pyogenes, Salmonella choleraesuis, Listeria monocytogenes,
	EN13727	1 minute	MRSA, Klebsiella pneumonia, Staphylococcus aureus, Pseudomonas aeruginosa, Enterococcus hirae, Acinetobacter baumannii
	EN13697	1 minute	Pseudomonas aeruginosa, Staphylococcus aureus, Escherichia coli, Enterococcus hirae
	EN1500*	30 seconds	Escherichia coli
Fungicidal	EN1650	1 minute	Candida albicans, Aspergillus niger, Tricophyton mentagrophytes
	EN13624	1 minute	Aspergillus niger, Candida albicans
Tuberculocidal	EN14348	1 minute	Mycobacterium terrae, Mycobacterium avium
Virucidal	EN14476	1 minute	HIV, Human Corona virus, Human Influenza virus, Herpes simplex virus, Swine Influenza A virus
	DVV/RKI**	1 minute	Vaccinia virus, BVDV (Hepatitis C surrogate)
Sporicidal	EN13704	1 minute	Clostridium difficile spores, Bacillus subtilis spores, Clostridium sporogenes spores, Bacillus cereus spores

\* The EN1500 test protocol (Chemical disinfectants and antiseptics. Hygienic handrub) was carried out on the hand wipe and hand foam presentations.

\*\*According to the DVV & RKI testing criteria, Spore-EX can be classified as active against all enveloped viruses in a 2% concentration in 1 minute.

### PromoCell GmbH

Sickingenstr. 63/65  
69126 Heidelberg  
Germany

Email: [info@promokine.info](mailto:info@promokine.info)

[www.promokine.info](http://www.promokine.info)

### North America

Phone: 1 – 866 – 251 – 2860 (toll free)  
Fax: 1 – 866 – 827 – 9219 (toll free)

### Deutschland

Telefon: 0800 – 776 66 23 (gebührenfrei)  
Fax: 0800 – 100 83 06 (gebührenfrei)

### France

Téléphone: 0800 90 93 32 (ligne verte)  
Téléfax: 0800 90 27 36 (ligne verte)

### United Kingdom

Phone: 0800 – 96 03 33 (toll free)  
Fax: 0800 – 169 85 54 (toll free)

### Other Countries

Phone: +49 6221 – 649 34 0  
Fax: +49 6221 – 649 34 40